Basic knowledge for catheter-based intervention

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THE CATHOLIC UNIVERSITY OF KOREA
INCHEON ST. MARY'S HOSPITAL

Angiography Systems

- **Fixed angiography system**
- **Mobile C-arm system**
- **By Special Function / Purpose**
 - **Hybrid OR angiography system**
 - **3D rotational angiography (3D-RA) system**
 - **Cone-beam CT angiography system**
 - **Fusion imaging / roadmap system**

Fluoroscopy Systems

C arm system

Hybrid system




Fig. 1.7.1.1 Diagram shows the components of a fluoroscopic imaging chain.

RADIOGRAPHIC IMAGING

- **Image Intensifier : 15"**
- **Digital Subtraction Angiography vs. Cineangiography**
 - Only the contrast-filled vascular structures(DSA)
 - Multiple x-ray pictures of the contrast-filled vessel as well as the surrounding tissue(Cine)
- **Trace Subtract Fluoroscopy or Road Mapping**
 - DSA +fluoroscopic overlay image
- **Leg Imaging Capability**
 - lower extremity angiography : head to feet
 - The table should be sufficiently long
 - Long wires (260 cm to 300 cm) are frequently used

Cineangiography vs. Digital Subtraction Angiography

Road map(Trace subtract fluoroscopic image)

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


Table 19-3

Recommended Radiographic Filming Projections for Optimal Branch Separation

Location	Recommended Filming Projection
Aortic arch	30-degree left anterior oblique
Cervical carotids	AP, lateral, and 45-degree ipsilateral anterior oblique
Intracranial carotids	AP and lateral
Vertebrobasilar system	AP and lateral
Right subclavian	Right anterior oblique
Renal artery origins	AP \pm 10 degrees
Celiac artery and SMA origins	Lateral
Iliac bifurcation	20- to 30-degree contralateral anterior oblique
Femoral bifurcation	20- to 30-degree ipsilateral anterior oblique
Trifurcation and tibial arteries	Anatomic AP (or 20-degree ipsilateral anterior oblique with feet in the neutral supine position)



AP, Anteroposterior; SMA, superior mesenteric artery.

Aortic arch tangential angle(θ -AATA)



Proximal Landing Zone

Position the image intensifier **LAO** for optimal imaging of aortic arch proximal landing zone

Distal Landing Zone

Position the image intensifier in **steep lateral** for optimal imaging of distal landing zone

Contrast agents (Iodinated contrast media)

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- CKD 환자: 3-4 mL/kg

Devices for Injection of Contrast Agents

- Power injection
- Hand injection

Table 19-2		Typical Injection Method, Rate, and Volume of Contrast Agent for Various Vascular Regions	
Location	Suggested Method	Injection Rate (mL/s)	Total Volume (mL)
Aortic arch	Power injection	20	40
Selective carotid	Hand or power injection	3-5	5-10
Selective vertebral	Hand injection	2-4	2-4
Selective subclavian or brachial	Hand or power injection	5-10	10
Abdominal aorta	Power injection	20	40
Renal or mesenteric	Hand injection	3-5	5-10
Iliac artery	Hand or power injection	10	10
Infrainguinal segments	Hand or power injection	5-10	10
Aorta to pedals, stepped run	Power injection	20	90*


Angiographic injector


ACCESS

- **One of the first steps, sono or fluoroscopic guided, micro-puncture**
- **Common Femoral Artery**
 - most common access site
 - sizes of 12 Fr to 14 Fr with minimal risk for ischemia
 - Most equipment has been developed for the femoral approach
 - Antegrade (toward the foot) or Retrograde (toward the iliac artery)

Optimal femoral artery access site

Micropuncture Set

Ideal artery puncture technique

- **Anatomy**
- **Fluoroscopic guided**
- **Sonoguided puncture**
- **Using micro puncture kit**
- **Optimal hemostasis method**

Arterial access

WIRE SELECTION (0.035", 0.018", 0.014")

■ 0.035"

- Subclavian, innominate, iliac, superficial femoral, and popliteal artery interventions.

■ 0.018 "

- Renal, subclavian, and popliteal interventions.

■ 0.014"

- Carotid, vertebral, renal, and below-knee interventions, coronary

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Wire Diameter	Directional Control	Tracking Capability	Support	Length (cm)	Company
<i>WIRES 0.035"</i>					
Stiff Angled Glidewire (3 cm or 5 cm flexible tip)	Excellent	Excellent	Good	260	Boston Scientific
Amplatz Extra Stiff w/ J-tip	Poor	Fair	Very Good	260	Boston Scientific
Amplatz Super Stiff (1 cm, 3.5 cm, 6 cm, and J-tip)	Poor	Poor	Excellent	260	Boston Scientific
Magic Torque	Fair	Good	Good	300	Boston Scientific
Wholey	Excellent	Good	Fair	260	Mallinckrodt
Lunderquist	Fair	Poor	Excellent	260	Cook
<i>WIRES 0.018"</i>					
Flex-T	Good	Good	Good	295	Mallinckrodt
Glidewire Gold (45 or 70 degree angle) with GT Leggiero Hydrophilic Microcatheter	Good	Excellent	Poor*	180	Boston Scientific
<i>WIRES 0.014"</i>					
Balance Trek	Good	Good	Fair	190, 300	Guidant

Table 1. Overview of different guidewires (0.018" and 0.014") and characteristics in endovascular peripheral arterial interventions

Manufacturer	Workhorse wire Stenotic/fibrotic lesions	CTO wire Calcified lesions	CTO wire Heavily calcified lesions
Asahi*	<i>Gladius 0.014/0.018</i> Sliding wire (intraluminal) <i>Gladius MG 0.014/0.018</i> Dissecting wire with microgap (subintimal)	<i>Halberd 0.014/0.018 (12 gram)</i> Drilling calcified lesion (preshape option) <i>Gaia PV 18 (7.5 gram)</i> Drilling lesion intraluminal	<i>Astato 30 0.018 (30 gram)</i> <i>Astato 40 0.014 (40 gram)</i> Tapered tip Penetrating and drilling calcified lesion
Abbott†	<i>Command 0.014/0.018 ST</i> Fibrotic lesion/looping (subintimal) <i>Command ES 0.014</i> More supportive variant (subintimal)	<i>Proceed 220T 0.014 (14 gram)</i> Preshape drilling calcified lesion <i>Winn 200T 0.014 (14.3 gram)</i> Tapered tip/Penetrating calcified lesion	<i>Connect 250T 0.018 (38.1 gram)</i> Tapered tip Penetrating calcified lesion
Boston Scientific‡	<i>V14/V18 ST</i> Loop/knuckle wire (subintimal)	<i>Victory 14 (12 gram)</i> <i>Victory 18 (12 gram)</i> penetrating calcified lesion	<i>Victory 14 (30 gram)</i> <i>Victory 18 (30 gram)</i> Penetrating calcified lesion
Terumo§	<i>Glidewire Advantage 14/18</i> Loop/knuckle wire (subintimal)	<i>Glidewire Advantage Track 14/18</i> Loop/knuckle wire (subintimal)	<i>Glidewire Advantage Track 14/18</i> Loop/knuckle wire (subintimal)
Cook	<i>Roadrunner Uniglide 0.018</i> <i>Approach Hydro ST 0.014</i> Loop/knuckle wire (subintimal)	<i>Approach CTO 0.014 (12 gram)</i> Penetrating calcified lesion	<i>Approach CTO 0.014 (24 gram)</i> Penetrating calcified lesion

CTO = chronic total occlusion; PV = peripheral vessel; ST = short tapering; ES = extra support.

* Asahi Intecc (Seto-shi, Aichi, Japan).


† Abbott (Chicago, IL, USA).

‡ Boston Scientific (Marlborough, MA, USA).

§ Terumo (Shibuya City, Tokyo, Japan).

|| Cook (Bloomington, IN, USA).

Wires...

Lunderquist Wire Guides 201Cm이상 (William Cook Europe)	0.35"*260 Cm	1	EA
Micropuncture Push-Plus Introducer Set Nitinol 4.0Fr*10Cm MPIS-401-NT-SST (Cook Medical)	Nitinol 4.0Fr*10Cm	1	EA
Performa Catheter All Size (Merit Medical)	JR*4.0*5Fr	1	EA
Radifocus Guide Wire 200Cm미만 (Terumo)	035*150*3 *A	1	EA
Radifocus Guide Wire 200Cm이상 (Terumo)	035*260*3 *A	1	EA
Radifocus Introducer II All Size, 4Fr~8Fr/7Cm, 10Cm (Terumo Corporation)	5Fr, 10Cm, A type	1	EA
Radifocus Introducer II All Size, 4Fr~8Fr/7Cm, 10Cm (Terumo Corporation)	7Fr, 10Cm, A type	1	EA

DIAGNOSTIC CATHETERS

■ Side-hole Diagnostic Catheters

- Permit large volumes of contrast agent to be infused safely in a large artery at a rapid rate via power injection
- Omniflush, Grollman, Universal Flush, Multipurpose, Pig tail

■ End-hole Diagnostic Catheters

- Selective diagnostic angiography is performed using manual injections of contrast agent through end-hole catheters
- 5Fr JR catheters...

Sample of side-hole catheters

A.Omni Flush


B.Multipurpose


C.Pigtail

Guide catheters

A: Headhunter (H-1); **B:** Judkins Right (JR4); **C:** Renal Standard (RES); **D:** Renal Multipurpose; **E:** Renal Double Curve (RDC); and **F:** Amplatz Right (AR1).


End-hold catheters

A: Cobra;
B: IMA;
C: Judkins Right
D: Simmons;
E: SoS
F: angled glide;
G: straight glide;
H: Vitek;
I: Berenstein;
J: red rubber.

Catheters...

Aortic advance catheter

David, Jerkin catheter



Angiographic catheter

Pigtail (5Fr)



SHEATH

- Using either a sheath or guide-based system
- Straight or shaped sheaths are most commonly used.
- These serve as the access points for delivery of equipment
- A hemostatic valve allows blood to bleed back when the valve is open, thereby reducing the risk of air and atherosclerotic debris embolization.

Introducer sheath

Sheaths...



A: and **B:** Low and high magnification view of Ansel sheath; **C:** and **D:** Low and high magnification view of Raabe sheath; **E:** and **F:** Low and high magnification view of Renal Double Curve guide sheath; and **G:** and **H:** Low and high magnification view of Shuttle sheath

Telescoping Technique

- **Advancing large guiding catheters and sheaths into certain arteries** may often be quite challenging and is often made simpler by telescoping the larger guide or sheath over a smaller diameter catheter.
- It is **quite important to have a smooth transition** between the guide or sheath and the telescoping catheter when advancing these devices around corners, and in tortuous diseased vessels, in order **to minimize vessel trauma and embolization**

Telescope technique


BALLOONS

- Most peripheral balloons are compliant with rated-burst pressures somewhat lower than their coronary counterparts
- To help these balloons fill and empty more effectively, a mixture of 70% saline and 30% contrast agent is used.
- The longer balloons are important when performing an SFA intervention from the brachial access site.

Balloon and inflator



Balloon의 역할


- 병변 확장으로 혈류 회복
- Stent 삽입 전 병변 준비 (pre-dilatation)
- Stent 삽입 후 확장 (post-dilatation)
- Drug delivery (Drug-Coated Balloon, DCB)

Balloon의 분류

- POBA (Plain Old Balloon Angioplasty)
- Drug-Coated Balloon (DCB)
- Cutting / Scoring Balloon
- High-pressure / Non-compliant Balloon
- Special Balloons (Long, Short)

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STENTS: BALLOON EXPANDABLE VS. SELF-EXPANDING, STENT GRAFTS

■ Balloon-expandable stents

- Require positive pressure for expansion and are typically rigid with high radial strength
- **Ideal for immobile areas of the body**, such as the subclavian, renal, mesenteric, and iliac arteries, and at ostial locations.

■ Self expandable stent (flexible stent)

- Nitinol stents : **flexible and memory** , Pulling back the sheath to allow stent expansion
- Carotid, axillary, and superficial femoral flex and twist whenever the neck, shoulder, and leg move

■ Stent graft (covered stent graft)

- **Exclude aneurysms and to treat perforations** when other more conservative options fail

Balloon Expandable Stent vs. Self-Expandable Stent

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
Vascular Closure Devices

- Adequate hemostasis after sheath removal is critical to avoid vascular and bleeding complications
- Manual compression : smaller arterial sheaths
- Suture-Based VCD : Perclose Proglide
- Collagen Plug : Angio seal, femo seal...
- Clip-Based
- Sealant-Based

Closing device

Thank you for your attention

THE CATHOLIC UNIVERSITY OF KOREA
INCHEON ST. MARY'S HOSPITAL